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## Instruction Manual

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### Material-related characteristics

#### Chemical Composition [wt. %]

Cobalt (Co)	63
Chromium (Cr)	29
Molybdenum (Mo)	6
Mn, Nb, Si, Fe	< 1

#### Physical specifications

Density	[g/cm <sup>3</sup> ]	8.8
E modulus	[GPa]	250
CTE (77–932°F)	[10 <sup>-6</sup> K <sup>-1</sup> ]	14.4
0.2 % proof stress	[MPa]	400
Vickers hardness	[HV10]	285
Elongation at break	[%]	9
Tensile strength	[MPa]	680
Melting point/- range:	[°C]	1350 - 1390

### Explanation of symbols

Manufacturer

Date of manufacture

Consult instructions for use

Use by

Batch

Catalogue No.

Store dry

Protect from direct sun light

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

#### 1. Indications

DD Bio CW are CoCr milling blanks for the manufacture of inlays/onlays, crowns & bridges of any span, primary & secondary parts, cone & telescopic crowns and bars & implant-supported superstructures.

#### 2. Contraindications

For patients who bear a potential risk of allergic reactions to one or more of the metals contained in the alloy, carrying out an allergy test prior to the dental treatment to guarantee the tolerance of the product is recommended.

#### 3. Safety information

Please pay attention to the information in the safety data sheet in its current version. Always wear suitable protective clothing while machining metals, and activate suction when working.

#### 4. Storage conditions

There are no special storage conditions.

#### 5. Preparation

For processing DD Bio CW Blanks, it is required to use appropriate metal milling tools and set the cooling lubricant at approx. 15%.

Please use at blanks with a high amount 15mm cutter with a larger exemption.

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### 6. Separating the work from the blank

To separate the frames from the blank, the use of suitable carbide cutters or appropriate cutting wheels is recommended.

### 7. Finishing / Cleaning

Finish the separated frames according to their intended application using carbide cutters or diamond burs. The material removal process should be carried out in such a way that the tools are only moved in one direction and material overlaps are prevented.

Observe the manufacturers' recommended maximum milling tool speed during the finishing process.

For your own safety, avoid the generation of excessive heat during the metal machining process.

After the conclusion of the working process with the hand piece, the frame surfaces should be grit-blasted with aluminium shot and the entire frame cleaned with steam or in distilled water in an ultrasonic device. A subsequent degreasing process with ethyl alcohol is recommended.

After cleaning, hold the frame only using an artery clamp or other appropriate tool.

### 8. Connecting

To connect the elements made of DD Bio CW, a cobalt-chrome solder or conventional cobalt-chrome laser welding wire should be used.

### 9. Oxidising

The implementation of oxidation firing is not necessary. If this should nevertheless be desired we recommend a temperature of 1000 °C (1832 °F) under atmospheric pressure and an immediate cooling is recommended. A further blasting process (aluminium oxide; approx. 180 µm; 3 - 4 bar) can afterwards be carried out. After this the surface has to be cleaned again.

### 10. Firing the ceramic

Conventional ceramic materials for dental alloys with a suitable CTE (< 14.4 x 10<sup>-6</sup>/K (25-500°C)) and veneering resin are suitable for veneering.

Please observe the working instructions from the respective manufacturer, in particular the cooling conditions specifically in massive constructions.

### 11. Polishing

First, remove the residual oxides with a polishing blasting agent and smooth out the coarse metal surfaces with a ceramic-bonded grindstone. Then polish all the relevant surfaces with a rubber polisher until they have a silk-gloss sheen. Finally, polish the metal surfaces with a polishing paste until they have a high-gloss sheen and then place the dental prosthesis in the ultrasonic device or clean it carefully with a steam-jet device.

### 12. Cementing

For temporary cementation are provisional cements or glass ionomer cements, recommended for permanent self-adhesive cementation cements.

### 13. Guarantee

All of the published recommendations from Dental Direkt GmbH are based on our own experience. The user is responsible for the proper application and working process.

Please find the current version of the instruction manual on our website: [www.dentaldirekt.de](http://www.dentaldirekt.de)

**This version replaces all previous versions.**